



# ENT/Neurosurgery

## High and low speed

### Rotary instruments

Komet Medical offers rotary instruments for many different hand-pieces and power systems.

Our range includes:

- Round burs made of stainless steel and tungsten carbide
- Diamond abrasives
- Cylindrical abrasives for special applications such as cochlear implantations

High-performance craniotomes for use in neurosurgery are also included in our range. To facilitate their use, our catalog and packaging label specify exactly which attachments our instruments are compatible with. Again, sterile packed instruments are available on request. Our rotary instruments are normally reusable and delivered a non-sterile condition.



